


“Unravel”


An ImageJ Macro to aid conversion of interface roughness to planar roughness

Peter J. Lee

Applied Superconductivity Center
NHMFL, Florida State University



Issue: How do you obtain roughness from a cross-section of an interface?





- Image above is a filament cross-section with a rough interface
 - How can the surface roughness be measured?



- Image above is an interface
 - Can the interfacial roughness be measured?



Menu 1: Choose object type

- Current version works on continuous outlines or solid objects
- The macro guesses the object type from measured parameters
 - In this case a horizontal line




Menu 2 (line): Sample length and basic options

- 3 options for providing measuring the total evaluation length are provided:
 - This evaluation length will be embedded in the exported map
- The outline can be skeletonized to minimum pixel width but at the risk of losing significant pixels. Fewer skeletonized pixels will be missed from the unraveled sequence, so it is the default setting.
- A pseudo height map can be exported that can be used in 3rd party topographic software like Gwyddion



Menu 3 (line): Unravelling options

- The starting point of the unravelling can be selected
- The search range for the next pixel can be adjusted
- A map can be created to show the pixel sequence order and the spline fit


A screenshot of a software dialog box titled "Unraveling options 3: (Unravel_interface_v221108.ijm)". The dialog has a light gray background and a dark gray title bar. It contains the following elements:

- Starting point:** A radio button is selected next to "Left pixel (4,55)". Below it is an unselected radio button for "Manual entry".
- Manual start x:** A text box containing the number "0" followed by the unit "pixels".
- Manual start y:** A text box containing the number "0" followed by the unit "pixels".
- Pixel search range in plus and minus pixels:** A text box containing the number "6" followed by the unit "pixels".
- Create map:** A checked checkbox with the label "Create map to show pixel sequence and spline fit".
- Buttons:** "OK" and "Cancel" buttons at the bottom right.



“Heights” from horizontal or vertical lines

- For horizontal or vertical lines, “heights” are simply obtained from the y coordinates or x coordinates respectively
- For horizontal lines the y coordinate is subtracted from the maximum y so that increasing “height” is upwards




	Seq_coord_x	Seq_coord_y	Seq_dist(px)	Seq_dist(μm)	Height(px)	Height_norm(px)	Height_norm(μm)	Height_norm^2(μm^2)	Fourier_Ht
1	4	65	0.000000000	0.000000000	33	13	5.038761701	25.389119484	16.099252701
2	4	64	1.000000000	0.387597054	34	14	5.426358755	29.445369342	11.786406517
3	3	64	2.000000000	0.775194108	34	14	5.426358755	29.445369342	1.603672147
4	5	64	4.000000000	1.550388216	34	14	5.426358755	29.445369342	1.291527987
5	6	64	5.000000000	1.937985270	34	14	5.426358755	29.445369342	0.290924072

Note:

The Fourier sequence is calculated for the entire sequence (column is for convenience).


Ra and Rq are not adjusted for true sampling lengths



File	Edit	Font
For 158B_67a_acBleach_(red)-skel_lzw.tif: Sub-sample spline-fit = 182.6171 μm		
158B_67a_acBleach_(red)-skel height statistics: min = 0 μm max = 26.3566 μm range = 26.3566 μm mean = 16.4094 μm std Dev = 6.3754 μm Ra(full length) = 16.4094 μm Rq(full length) = 17.6029 μm		


Pseudo height map from horizontal or vertical lines

- The 2D line data can be stretched vertically to create a 32-bit pseudo height map that can be imported into 3rd party topological analysis software
- Unit-scaled “Heights” are directly stored in the 32-bit data
- The measured evaluation length is embedded in the horizontal scale





Additional output

- If the sequence map option is chosen the spline fit is drawn over the map as an overlay.





Menu 1: Choose object type – outline continuous

- This example: A continuous outline around an object





Menu 2 (continuous): Sample length and basic options

- 5 options for providing measuring the total evaluation length are provided:
- This evaluation length will be embedded in the exported map



Menu 3 (continuous): Unravelling options

- The starting point of the unravelling can be selected
- The search range for the next pixel can be adjusted
- A map can be created to show the pixel sequence order and the spline fit



Unraveling options 3: (Unravel_interface_v221110.ijm) X

Starting point:

☒ Left pixel (123,657)

☐ Top pixel (615,48)

☐ Manual entry

Manual start x pixels

Manual start y pixels

Pixel search range in plus and minus pixels pixels

☒ Try to start clockwise?

☒ Output rotational sequence values

☒ Create map to show pixel sequence and spline fit

OK Cancel

“Height” reference menu (continuous)

- For continuous objects, distances (“heights”) are calculated to a single reference point (default is the measured object center)

Height reference coordinate (Unravel_interface_v221110.ijm)

Sequential pixel centroid: x = 594.0077, y = 545.4681
Object center: x = 589.5622, y = 515.9225
Image center: x = 640, y = 512
Reference shape (Ellipse) center: x = 566.5, y = 524

Reference location for height:

☐ Sequential_pixel_centroid
☒ Object_center
☐ Image_center
☐ Reference_shape_center
☐ Arbitrary_coordinates

Arbitrary x: pixels
Arbitrary y: pixels

Table compatible name for distance to reference (i.e. 'Height'):


OK Cancel

- The minimum distance value is subtracted to output the normalized values

Results												
	Seq_coord_x	Seq_coord_y	Seq_dist(px)	Seq_dist(μm)	Height(px)	Height_norm(px)	Height_norm(μm)	Height_norm^2(μm^2)	Angle (radians)	Angle Offset (radians)	Angle Offset (degrees)	Fourier_amps
1	123	657	0.000000000	0.000000000	487.424964538	129.981082914	50.380303430	2538.174973736	0.185072240	0.000000000	0.000000000	28.321575165
2	123	656	8.544003745	3.311631905	487.136470741	129.692589116	50.268484043	2526.920487976	0.185347950	0.000275710	0.015797016	11.625046730
3	123	655	9.544003745	3.699229102	486.849860014	129.405978390	50.157394529	2515.764225915	0.185624473	0.000552233	0.031640627	13.900227547
4	123	654	10.544003745	4.086826299	486.565135687	129.121254063	50.047036178	2504.705830165	0.185901813	0.000829573	0.047531037	3.644485712
5	124	653	11.958217308	4.634971512	485.322935813	127.879054188	49.565562988	2456.745034297	0.187658482	0.002586242	0.148180744	6.416252613

Additional output

- A basic summary is sent to the log window
- Note: The Fourier amplitudes are calculated for the entire sequence (column is for convenience).
- Ra and Rq are not adjusted for true sampling lengths



```
Log
File Edit Font

For 158B_67a_acBleach(red)_OtsuThreshInLine_lzw.tif:
Ellipse perimeter = 1050.6251  $\mu\text{m}$ 
Reference locations: x = 589.5622, y = 515.9225

158B_67a_acBleach(red)_OtsuThreshInLine_unrav height statistics:
min = 0  $\mu\text{m}$ 
max = 61.161  $\mu\text{m}$ 
range = 61.161  $\mu\text{m}$ 
mean = 28.3216  $\mu\text{m}$ 
std Dev = 10.7353  $\mu\text{m}$ 
Ra(full length) = 28.3216  $\mu\text{m}$ 
Rq(full length) = 30.2876  $\mu\text{m}$ 

Warning: 1060 non-Sequenced pixels
```


Note: If not all pixels are captured in the sequence the missing count as reported and an image (left) is created with the missing pixels highlighted. The pixel capture is not uniform because the search matrix looks for the first-nearest pixel with the search sequence the same for all searches

Some improvement might be useful here

The higher the image resolution for the interface, the less likelihood of error from this issue



Summary Image

- A summary image is created showing:
 - Direction of sequence (dark-to-light)
 - Fit used for evaluation length (green)
 - Any missed pixels from the sequence (surrounded by blue highlight)






Pseudo height map from horizontal or vertical lines

- The 2D line data can be stretched vertically to create a 32-bit pseudo height map that can be imported into 3rd party topological analysis software.
- Unit-scaled “Heights” are directly stored in the 32-bit data
- The measured sample length is embedded in the horizontal scale



Gwyddion Import (<http://gwyddion.net/>)

- Drop TIFF format hMap file onto Gwyddion toolbar to open import window




Leave as “1” to directly import 32-bit height data from the TIFF file


The last width is remembered so it is important to correct this to the current value. Should this value be the circle/ellipse perimeter or the smoothed surface?

These are the units used for the 32-bit TIFF file

Import into Gwyddion to correct curvature and measure roughness parameters



Roughness values can be calculated by clicking on this icon on the Gwyddion toolbar



- After clicking on icon draw a horizontal line across the map (use shift key to enforce horizontal)

